



heavy duty nursing home bed Gigant



nursing home bed Gigant

> performance

4-section mattress base 120 x 200 cm, mattress surface height adjustment from 40 to 80 cm, back rest and thigh rest as well as the height of the mattress surface are adjustable by electric motors. The foot section height of the mattress base can be lowered by tilting the mattress base (adjustment by hand control). Metal parts are sealed with a polyester powder lacquer coating, in the colour topaz. The wooden surround is in beech finish/solid beech.

safe working load: 350 kg
max. patient's weight: 280 kg

manufactured according to:

- > norm DIN EN 1970:2000 + A1:2005
- > norm DIN EN 60601-2-38/A1:2001 (reduced for home care requirements)
- > norm IEC 60601-2-52:2009

> order numbers

- | | |
|-----------|--|
| 59.6631.1 | nursing home bed "Gigant"
(wooden side rails) |
| 19.1400.1 | lifting pole
(included) |

> dimensions

length 220 cm, width 134 cm

> supplementary equipment (available at extra cost):

(technical changes reserved)



clamp-on extra-high side rail



reading lamps




integrated bed extension including 4 longer wooden side rails 220 cm (also for retrofitting)



padded side rail cover



Gigant - the nursing home bed for heavyweight residents up to 280 kg body weight



350 kg safe working load, 4-section mattress base, electrically adjustable

drive system DEWERT being easy to maintain and durable

mattress base extremely robust with steel grating floor, mattress base width: 120 cm

extremely solid construction of side rails

locking device included

clearance with all commercially available lifters



lowered foot section



lower leg rest liftable by ratchet

your advantages:

- **quality and durability**

The home care bed "Gigant" is extremely solid constructed, so that it fulfills the requirements in daily care perfectly. The high quality ensures you a nursing home bed which is durable, robust and easy to maintain and which decisively minimizes the follow-up costs.

- **350 kg safe working load**

The mattress base of the nursing home bed "Gigant" can accommodate a safe working load up to 350 kg (body weight + mattress + bedding + lifting pole etc. < 350kg). Residents weighing up to 280 kg find the best possible conditions on the 120 cm wide mattress base. If necessary, the mattress base can be extended to 220 cm (bed extension available as supplementary equipment).